


High-Throughput Measurement of Deep Eutectic Solvent Melting Points using IR Bolometry


Project Sponsor: Dr. Lilo D. Pozzo (*Chemical Engineering*)

Team Members: Shrilakshmi Bonageri (*Chemical Engineering*), Jaime Rodriguez (*Chemical Engineering*), Sage Scheiwiller (*Chemical Engineering*)




Background and Motivation


Deep eutectic solvents (DES) are novel solvents formed between organic hydrogen bond acceptors and donors. DES can be formed at low-cost for several important applications, such as chemical synthesis, extractions, electrochemistry, and even drug delivery. However, the design space for DES is enormous and high-throughput measurement of melting points is required to rapidly identify DES with melting points that are feasible for their intended application.




Workflow






Bolometry Setup




- Sample temperature is monitored by an IR camera.
- Melting points are detectable from a sudden increase in sample temperature due to an increased thermal conductivity in the liquid phase.



Temperature Profile via Edge Detection




- The temperature profile of the samples and plate is determined by detecting the edges, filling and labeling them, and monitoring the temperature at their centroids.




Alternative Method for Low-Contrast Samples

- In some situations, the contrast between the image and sample may be too low for edge detection, even with contrast enhancement.





- Alternatively, centroid locations for each sample can be found by summing pixel values over individual rows and columns of the sample holder (well plate).

Sum of pixel values over columns against column index



Sum of pixel values over rows against row index



- Using this alternative method, centroids were located for each sample in a 96 well microplate.

Conclusions and Future Work

Using IR Bolometry, melting points can be determined for multiple samples at once in a matter of minutes, as opposed to standard techniques which may take up to an hour for a single sample. Once sufficient data is collected, future work may include the development of a machine learning model to predict the melting points of DES based on their composition.